


Anti-Tubgcp6 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Rt7

Item Number: CM0018976



Introduction

Rabbit polyclonal antibody targeting TUBGCP6 (GCP6), a key gamma-tubulin ring complex subunit. Suitable for WB, IHC-P, ELISA in human, mouse, rat. Essential for centrosome and microtubule nucleation studies.

[Learn More](#)

Parameter	Value
Product Name	TUBGCP6 Rabbit pAb
Remarks/Alias	GCP6; GCP-6; MCCRP; MCCRP1; MCPHCR; TUBGCP6
Species	Human
Gene ID (Human)	85378
Gene ID	85378
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-900 of human TUBGCP6 (NP_065194.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96RT7
Protein Weight	200kDa
Shipping	Ice bag