

Anti-Tubulin Beta-1 Chain Monoclonal Antibody For Wb, If/Icc, Elisa - Q9H4B7

Item Number: CM0002060

Introduction

Anti-Tubulin beta-1 chain rabbit monoclonal antibody for WB, IF/ICC, ELISA. Recognizes human Tubulin beta-1 chain, cross-reacts with mouse, rat, cow, dog, and monkey. Suitable for microtubule structure and function studies.

[Learn More](#)

Parameter	Value
Product Name	Tubulin beta-1 chain Rabbit mAb
Remarks/Alias	MACTHC1; Tubulin beta-1 chain
Species	Human
Gene ID (Human)	81027
Gene ID	81027
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 352-451 of human Tubulin beta-1 chain (Q9H4B7).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Cow, Dog, Monkey
SWISS	Q9H4B7
Protein Weight	50 kDa
Shipping	Ice bag