


Anti-Tufm Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P49411

Item Number: CM0025723



Introduction

High-quality rabbit polyclonal antibody against human TUFM (EF-TuMT), validated for WB, IHC-P, IF/ICC, IP, ELISA. Ideal for mitochondrial translation, autophagy, and innate immunity research. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	TUFM Rabbit pAb
Remarks/Alias	P43; EFTU; COXPD4; EF-TuMT; TUFM
Species	Human
GeneID (Human)	7284
GeneID	7284
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 186-455 of human TUFM (NP_003312.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49411
Protein Weight	50kDa
Shipping	Ice bag