


Anti-Tug/Aspscr1 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Bze9

Item Number: CM0003903



Introduction

CamelBio's Anti-TUG/ASPSCR1 rabbit monoclonal antibody, validated for Western Blot and ELISA. Detects human and mouse TUG/ASPSCR1, a GLUT4 tethering protein at ~60kDa. Ideal for metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TUG/ASPSCR1 Rabbit mAb
Remarks/Alias	TUG; ASPL; ASP5; RCC17; UBXD9; UBXN9; ASPCR1; TUG/ASPSCR1
Species	Human
Gene ID (Human)	79058
Gene ID	79058
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human TUG/ASPSCR1 (Q9BZE9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BZE9
Protein Weight	60kDa
Shipping	Ice bag