

Anti-Tulp2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P46686

Item Number: CM0020247

Introduction

Rabbit polyclonal antibody targeting mouse Tulp2 (Tubby-related protein 2). Validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Suitable for retina and testis research.

[Learn More](#)

Parameter	Value
Product Name	Tulp2 Rabbit pAb
Remarks/Alias	Pdet, Tulp2
Species	Mouse
Gene ID (Human)	7288
Gene ID	56734
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 78-305 of mouse Tulp2 (NP_032833.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46686
Protein Weight	63kDa
Shipping	Ice bag