

Anti-Tweakr/Fn14/Cd266 Polyclonal Antibody For Wb - Q9Np84

Item Number: CM0013660



Introduction

Rabbit polyclonal antibody targeting human TWEAKR/Fn14/CD266 (TNFRSF12A). Validated for Western blot and ELISA applications. Recombinant immunogen; 14 kDa; human reactive. Suitable for apoptosis and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	TWEAKR/Fn14/CD266 Rabbit pAb
Remarks/Alias	FN14; CD266; TWEAKR; TWEAKR/Fn14/CD266
Species	Human
GeneID (human)	51330
GeneID	51330
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-129 of human TWEAKR/Fn14/CD266 (NP_057723.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NP84
Protein Weight	14kDa
Shipping	Ice bag