

Anti-Tweakr/Fn14/Cd266 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Np84

Item Number: CM0010392



Introduction

Rabbit monoclonal antibody targeting human TWEAKR/Fn14/CD266 (TNFRSF12A), validated for WB and ELISA. Cross-reacts with mouse. Recombinant immunogen from amino acids 50-129. Recognizes 14kDa protein. Suitable for apoptosis, angiogenesis, and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	TWEAKR/Fn14/CD266 Rabbit mAb
Remarks/Alias	FN14; CD266; TWEAKR; TWEAKR/Fn14/CD266
Species	Human
GeneID (Human)	51330
GeneID	51330
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 50-129 of human TWEAKR/Fn14/CD266 (Q9NP84))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NP84
Protein Weight	14kDa
Shipping	Ice bag