

Anti-Twf1 Polyclonal Antibody For Wb, Elisa - Q12792

Item Number: CM0017720



Introduction

Rabbit polyclonal antibody targeting human TWF1 (Twinfilin-1/A6/PTK9). Validated for Western blot and ELISA. Cross-reacts with human and mouse. Ideal for actin dynamics and endocytosis research.

[Learn More](#)

Parameter	Specification
Product Name	TWF1 Rabbit pAb
Remarks/Alias	A6, PTK9, TWF1
Species	Human
Gene ID (Human)	5756
Gene ID	5756
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 101-350 of human TWF1 (NP_002813.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12792
Protein Weight	40kDa
Shipping	Ice bag