


Anti-Twf2 Polyclonal Antibody For Wb, Elisa - Q6lbs0

Item Number: CM0026428



Introduction

Rabbit polyclonal antibody against twinfilin-2 (TWF2), validated for Western blotting and ELISA. Reacts with human and mouse. Involved in actin dynamics and endocytosis. Suitable for cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	TWF2 Rabbit pAb
Remarks/Alias	A6r; A6RP; PTK9L; MSTP011; TWF2
Species	Human
Gene ID (Human)	11344
Gene ID	11344
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-349 of human TWF2 (NP_009215.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6lbs0
Protein Weight	40kDa
Shipping	Ice bag