


Anti-Twf2 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q6Ibs0

Item Number: CM0006961



Introduction

Rabbit monoclonal antibody targeting human TWF2 (Twinfilin-2), validated for WB, IP, and ELISA applications. Cross-reacts with human, mouse, and rat. Recognizes a ~40 kDa protein involved in actin dynamics and clathrin-mediated endocytosis.

[Learn More](#)

Parameter	Value
Product Name	TWF2 Rabbit mAb
Remarks/Alias	A6r; A6RP; PTK9L; MSTP011; TWF2
Species	Human
Gene ID (Human)	11344
Gene ID	11344
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-349 of human TWF2. (NP_009215.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6IBS0
Protein Weight	40 kDa
Shipping	Ice bag