

Anti-Twist Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15672

Item Number: CM0025904



Introduction

Rabbit polyclonal antibody against human Twist protein (TWIST1, bHLHa38). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for transcription regulation, craniofacial development, and myogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Twist Rabbit pAb
Remarks/Alias	CRS; CSO; SCS; ACS3; CRS1; BPES2; BPES3; SWCOS; TWIST; bHLHa38; Twist
Species	Human
GeneID (human)	7291
GeneID	7291
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of Twist (NP_000465.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15672
Protein Weight	21kDa
Shipping	Ice bag