

Anti-Twist Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15672

Item Number: CM0023779



Introduction

High-quality rabbit polyclonal antibody against human Twist (TWIST1) for WB, IHC-P, ELISA. Validated in human, mouse, rat. Ideal for studying transcription regulation, myogenesis, and craniosynostosis.

[Learn More](#)

Parameter	Value
Product Name	Twist Rabbit pAb
Remarks/Alias	CRS; CSO; SCS; ACS3; CRS1; BPES2; BPES3; SWCOS; TWIST; bHLHa38; Twist
Species	Human
Gene ID (Human)	7291
Gene ID	7291
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 51-150 of human Twist (NP_000465.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15672
Protein Weight	21kDa
Shipping	Ice bag