

Anti-Twist1 Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - Q15672

Item Number: CM0005905



Introduction

Validated TWIST1 rabbit monoclonal antibody for western blot, immunofluorescence, and ELISA. Reacts with human and mouse TWIST1, ideal for transcription regulation and developmental biology studies.

[Learn More](#)

Product Parameter	Value
Product Name	TWIST1 Rabbit mAb
Remarks/Alias	TWIST1; ACS3; BPES2; BPES3; CRS; CRS1; CSO; SCS; TWIST; bHLHa38; twist-related protein 1
Species	Human
GeneID (Human)	7291
GeneID	7291
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human TWIST1 (NP_000465.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15672
Protein Weight	21kDa
Shipping	Ice bag