

Anti-TxlNa Rabbit Polyclonal Antibody For Wb, Elisa - P40222

Item Number: CM0020154

Introduction

Anti-TXLNA rabbit polyclonal antibody validated for WB, ELISA. Recognizes human alpha-taxilin (62 kDa, gene ID 200081), a protein involved in intracellular vesicle trafficking and exocytosis. For research use.

[Learn More](#)

Parameter	Value
Product Name	TXLNA Rabbit pAb
Remarks/Alias	IL14; TXLN; TXLNA
Species	Human
Gene ID (Human)	200081
Immunogen	Recombinant protein corresponding to amino acids 1-160 of human TXLNA (NP_787048.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P40222
Protein Weight	62kDa
Shipping	Ice bag