

Anti-TxlNa Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P40222

Item Number: CM0019666

Introduction

Rabbit polyclonal antibody targeting human alpha-taxilin (TXLNA). Suitable for Western blot, IF/ICC, and ELISA. Ideal for studying intracellular vesicle traffic and calcium-dependent exocytosis.

[Learn More](#)

Parameter	Value
Product Name	TXLNA Rabbit pAb
Remarks/Alias	IL14; TXLN; TXLNA
Species	Human
GeneID (human)	200081
GeneID	200081
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human TXLNA (NP_787048.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P40222
Protein Weight	62kDa
Shipping	Ice bag