

Anti-Txndc12 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95881

Item Number: CM0018647



Introduction

Rabbit polyclonal antibody targeting human TXNDC12 (ERp18/ERp19), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying ER protein disulfide isomerase activity and redox regulation.

[Learn More](#)

Parameter	Value
Product Name	TXNDC12 Rabbit pAb
Remarks/Alias	AG1, AGR1, ERP16, ERP18, ERP19, TLP19, hAG-1, PDIA16, hTLP19, TXNDC12
Species	Human
GeneID (Human)	51060
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 93-172 of human TXNDC12 (NP_056997.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95881
Protein Weight	19kDa
Shipping	Ice bag