


Anti-Txndc5 Polyclonal Antibody For Wb - Q8Nbs9

Item Number: CM0018664



Introduction

Rabbit polyclonal antibody against human TXNDC5 (ERp46), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for studying ER protein folding and disulfide bond formation.

[Learn More](#)

Parameter	Value
Product Name	TXNDC5 Rabbit pAb
Remarks/Alias	HCC2; ERP46; HCC-2; STRF8; PDIA15; UNQ364; ENDOPDI; TXNDC5
Species	Human
GeneID (Human)	81567
GeneID	81567
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-324 of human TXNDC5 (NP_001139021.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NBS9
Protein Weight	48kDa
Shipping	Ice bag