

Anti-Txndc9 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O14530

Item Number: CM0030321

Introduction

Rabbit polyclonal antibody against TXNDC9, validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat TXNDC9, a negative regulator of chaperonin TCP1 ATPase activity involved in actin and tubulin folding.

[Learn More](#)

Parameter	Value
Product Name	TXNDC9 Rabbit pAb
Remarks/Alias	APACD, PHLP3, TXNDC9
Species	Human
Gene ID (Human)	10190
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human TXNDC9 (NP_005774.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14530
Protein Weight	27kDa
Shipping	Ice bag