


Anti-Txnip Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H3M7

Item Number: CM0026063



Introduction

CamelBio offers Rabbit anti-TXNIP polyclonal antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat TXNIP (~44kDa). Functions in oxidative stress, cell cycle arrest, and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	TXNIP Rabbit pAb
Remarks/Alias	THIF; VDUP1; ARDC6; HHCPA78; EST01027; TXNIP
Species	Human
Gene ID (Human)	10628
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-166 of human TXNIP (NP_006463.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H3M7
Protein Weight	44kDa
Shipping	Ice bag