


Anti-Txnl1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O43396

Item Number: CM0022924



Introduction

Rabbit polyclonal antibody targeting human TXNL1 (Thioredoxin-like protein 1), reactive with human, mouse, rat. Suitable for WB, IHC-P, IF/ICC, ELISA. Swiss-Prot O43396. Research use in redox biology.

[Learn More](#)

Parameter	Specification
Product Name	TXNL1 Rabbit pAb
Remarks/Alias	Txl; TXNL; TRP32; TXL-1; HEL-S-114; TXNL1
Species	Human
Gene ID (Human)	9352
Gene ID	9352
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-289 of human TXNL1 (NP_004777.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43396
Protein Weight	32kDa
Shipping	Ice bag