


Anti-Txnl4A Rabbit Polyclonal Antibody For Wb And Elisa - P83876

Item Number: CM0016977



Introduction

Polyclonal antibody targeting human TXNL4A (Thioredoxin-like protein 4A, also DIM1/U5-15kD). Validated for WB and ELISA with cross-reactivity to mouse. Suitable for spliceosome research.

[Learn More](#)

Parameter	Value
Product Name	TXNL4A Rabbit pAb
Remarks/Alias	BMKS, DIB1, DIM1, TXNL4, SNRNP15, U5-15kD, TXNL4A
Species	Human
GeneID (Human)	10907
GeneID	10907
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human TXNL4A (NP_006692.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P83876
Protein Weight	17kDa
Shipping	Ice bag