

Anti-Txnl4B Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nx01

Item Number: CM0028154



Introduction

Rabbit polyclonal antibody targeting TXNL4B (DLP/Dim2). Validated for WB, IF/ICC, ELISA. Cross-reacts with human and mouse. Involved in pre-mRNA splicing and S/G2 transition.

[Learn More](#)

Parameter	Value
Product Name	TXNL4B Rabbit pAb
Remarks/Alias	DLP; Dim2; TXNL4B
Species	Human
Gene ID (Human)	54957
Gene ID	54957
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-149 of human TXNL4B (NP_060323.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NX01
Protein Weight	17kDa
Shipping	Ice bag