


Anti-Tyk2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P29597

Item Number: CM0021533



Introduction

TYK2 Rabbit pAb detects human, mouse, rat TYK2 via Western blot, IHC-P, IF/ICC, and ELISA. This 134 kDa JAK family kinase mediates interferon and cytokine signaling, essential for innate and adaptive immunity.

[Learn More](#)

Parameter	Value
Product Name	TYK2 Rabbit pAb
Remarks/Alias	JTK1; IMD35; TYK2
Species	Human
Gene ID (Human)	7297
Gene ID	7297
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1102-1182 of human TYK2 (NP_003322.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29597
Protein Weight	134kDa
Shipping	Ice bag