

Anti-Tymp Polyclonal Antibody For Wb And Elisa - P19971

Item Number: CM0016134



Introduction

Rabbit polyclonal antibody targeting TYMP (Thymidine phosphorylase/PD-ECGF), validated for WB and ELISA in human and rat samples. Essential for angiogenesis and nucleoside metabolism research.

[Learn More](#)

Parameter	Value
Product Name	TYMP Rabbit pAb
Remarks/Alias	TP; ECGF; ECGF1; MNGIE; MEDPS1; MTDPS1; PDECGF; hPD-ECGF; TYMP
Species	Human
GeneID (Human)	1890
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P19971
Protein Weight	50kDa
Shipping	Ice bag