

Anti-Tyms Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P04818

Item Number: CM0018583



Introduction

TYMS Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Detects human and mouse thymidylate synthase (TS, TMS, HST422), a key enzyme in de novo thymidylate biosynthesis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TYMS Rabbit pAb
Remarks/Alias	TS; TMS; DKCD; HST422; TYMS
Species	Human
Gene ID (Human)	7298
Gene ID	7298
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-313 of human TYMS (NP_001062.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P04818
Protein Weight	36kDa
Shipping	Ice bag