

## Anti-Tyro3 Monoclonal Antibody For Wb And Elisa - Q06418

Item Number: CM0002124



### Introduction

Recombinant rabbit monoclonal antibody against human Tyro3 (Q06418), validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for studying cell survival, migration, and viral entry mechanisms.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | Tyro3 Rabbit mAb                                                                                                      |
| Remarks/Alias    | BYK; Dtk; RSE; Rek; Sky; Tif; Etk-2; Tyro3                                                                            |
| Species          | Human                                                                                                                 |
| Gene ID (Human)  | 7301                                                                                                                  |
| Gene ID          | 7301                                                                                                                  |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 791-890 of human Tyro3 (Q06418). |
| Source           | Rabbit                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                     |
| SWISS            | Q06418                                                                                                                |
| Protein Weight   | 97kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |