

Anti-Tyrosinase Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P14679

Item Number: CM0016599



Introduction

High-quality rabbit polyclonal antibody against human Tyrosinase, validated for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for melanin biosynthesis and pigmentation research.

[Learn More](#)

Parameter	Value
Product Name	Tyrosinase Rabbit pAb
Remarks/Alias	ATN; CMM8; OCA1; OCA1A; OCAIA; SHEP3; Tyrosinase
Species	Human
GeneID (Human)	7299
GeneID	7299
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 20-340 of human Tyrosinase (NP_000363.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14679
Protein Weight	60kDa
Shipping	Ice bag