

Anti-Tyrosinase Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P14679

Item Number: CM0018725



Introduction

Tyrosinase Rabbit Polyclonal Antibody (60kDa) suitable for WB, IHC-P, ELISA. Reacts with human, rat. Targets melanin biosynthesis enzyme; ideal for albinism and melanoma research.

[Learn More](#)

Parameter	Value
Product Name	Tyrosinase Rabbit pAb
Remarks/Alias	ATN; CMM8; OCA1; OCA1A; OCAIA; SHEP3; Tyrosinase
Species	Human
GeneID (Human)	7299
GeneID	7299
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-340 of human Tyrosinase (NP_000363.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P14679
Protein Weight	60kDa
Shipping	Ice bag