


Anti-Tyrosine Hydroxylase Monoclonal Antibody For Wb, If, Ihc-P, And Elisa - P07101

Item Number: CM0005365



Introduction

Mouse monoclonal antibody against human Tyrosine Hydroxylase (TH), suitable for WB, IF, IHC-P, and ELISA. Reacts with human, mouse, rat. Validated for research and diagnostic applications in catecholamine biosynthesis studies.

[Learn More](#)

Parameter	Value
Product Name	Tyrosine Hydroxylase Mouse mAb
Remarks/Alias	TYH; DYT14; DYT5b
Species	Human
GeneID	7054
Immunogen	A synthesized peptide derived from human Tyrosine Hydroxylase
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07101
Protein Weight	45-59 kDa
Shipping	Ice bag