


Anti-Tyrosine Hydroxylase Monoclonal Antibody For Wb, Ihc, If, Elisa - P07101

Item Number: CM0004093



Introduction

Recombinant rabbit monoclonal antibody targeting human Tyrosine Hydroxylase (TH, TYH). Validated for WB, IF/ICC, IHC-P, and ELISA, with cross-reactivity to mouse and rat. Suitable for catecholamine and Parkinson's disease research.

[Learn More](#)

Parameter	Value
Product Name	Tyrosine Hydroxylase Rabbit mAb
Remarks/Alias	TYH; DYT14; DYT5b
Species	Human
GeneID (Human)	7054
GeneID	7054
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 411-528 of human Tyrosine Hydroxylase (NP_954986.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07101
Protein Weight	59kDa
Shipping	Ice bag