


Anti-Tyrosine Hydroxylase Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P07101

Item Number: CM0012881



Introduction

High-quality Anti-Tyrosine Hydroxylase Rabbit Polyclonal Antibody validated for Western Blot, IF/ICC, ELISA in human, mouse, rat. Ideal for catecholamine biosynthesis research and IVD assay development.

[Learn More](#)

Parameter	Specification
Product Name	Tyrosine Hydroxylase Rabbit pAb
Remarks/Alias	TYH; DYT14; DYT5b; Tyrosine Hydroxylase
Species	Human
Gene ID (Human)	7054
Gene ID	7054
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 411-528 of human Tyrosine Hydroxylase (NP_954986.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07101
Protein Weight	59kDa
Shipping	Ice bag