


Anti-Tyrosine Hydroxylase Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P07101

Item Number: CM0016247



Introduction

Rabbit polyclonal antibody against human Tyrosine Hydroxylase, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets the rate-limiting enzyme in catecholamine biosynthesis, suitable for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Tyrosine Hydroxylase Rabbit pAb
Remarks/Alias	TYH; DYT14; DYT5b; Tyrosine Hydroxylase
Species	Human
GeneID (Human)	7054
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 411-528 of human Tyrosine Hydroxylase (NP_954986.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07101
Protein Weight	45-59 kDa
Shipping	Ice bag