

Anti-U1A/Snrpa Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09012

Item Number: CM0025341



Introduction

Rabbit polyclonal antibody targeting U1A/SNRPA, validated in WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for spliceosome and RNA processing research.

[Learn More](#)

Parameter	Value
Product Name	U1A/SNRPA Rabbit pAb
Remarks/Alias	U1A; Mud1; U1-A; U1A/SNRPA
Species	Human
GeneID (Human)	6626
GeneID	6626
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-282 of human U1A/U1A/SNRPA (NP_004587.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09012
Protein Weight	31kDa
Shipping	Ice bag