

Anti-U2Af1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q01081

Item Number: CM0017765



Introduction

Rabbit polyclonal antibody targeting human U2AF1, the 35 kDa subunit of the U2AF splicing factor. Proven in WB, IHC-P, IP, ELISA; cross-reacts with mouse and rat. Ideal for pre-mRNA splicing and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	U2AF1 Rabbit pAb
Remarks/Alias	RN; FP793; U2AF35; U2AFBP; RNU2AF1; U2AF1
Species	Human
Gene ID (Human)	7307
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human U2AF1 (NP_006749.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01081
Protein Weight	28kDa
Shipping	Ice bag