

Anti-U2AF2/U2AF65 Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P26368

Item Number: CM0009930



Introduction

Rabbit monoclonal antibody targeting human U2AF2/U2AF65 (P26368). Validated for WB, IHC-P, IP, ELISA; cross-reacts with human, mouse, rat. Ideal for splicing and mRNA processing studies.

[Learn More](#)

Parameter	Value
Product Name	U2AF2/U2AF65 Rabbit mAb
Remarks/Alias	U2AF65; U2AF2/U2AF65
Species	Human
GeneID (Human)	11338
GeneID	11338
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 376-475 of human U2AF2/U2AF2/U2AF65 (P26368).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26368
Protein Weight	54kDa
Shipping	Ice bag