

Anti-U2AF2/U2AF65 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P26368

Item Number: CM0022389



Introduction

Rabbit polyclonal antibody targeting human U2AF2/U2AF65. Validated for WB, IF/ICC, IP, and ELISA. Cross-reacts with human and mouse. Ideal for splicing, 3'-end processing, and mRNA export research.

[Learn More](#)

Parameter	Value
Product Name	U2AF2/U2AF65 Rabbit pAb
Remarks/Alias	U2AF65; U2AF2/U2AF65
Species	Human
GeneID (Human)	11338
GeneID	11338
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 192-471 of human U2AF2/U2AF65 (NP_001012496.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P26368
Protein Weight	54kDa
Shipping	Ice bag