

Anti-Uaca Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bzf9

Item Number: CM0028301

Introduction

High-quality rabbit polyclonal antibody against human UACA (NUCLING), validated for Western blot and ELISA. Reacts with human, mouse, and rat samples. Ideal for apoptosis and cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	UACA Rabbit pAb
Remarks/Alias	NUCLING; UACA
Species	Human
GenelD (Human)	55075
GenelD	55075
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZF9
Protein Weight	163kDa
Shipping	Ice bag