

Anti-Uap1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q16222

Item Number: CM0010899



Introduction

Rabbit monoclonal antibody targeting human UAP1 (Q16222), suitable for WB, IF/ICC, ELISA, with cross-reactivity to mouse and rat. Useful for investigating UDP-GlcNAc biosynthesis and innate immune signaling.

[Learn More](#)

Parameter	Value
Product Name	UAP1 Rabbit mAb
Remarks/Alias	AGX; AGX1; AGX2; SPAG2; UAP1
Species	Human
GeneID (Human)	6675
GeneID	6675
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 400-500 of human UAP1 (Q16222).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16222
Protein Weight	59kDa
Shipping	Ice bag