


Anti-Uap56/Ddx39B Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13838

Item Number: CM0007990



Introduction

Rabbit monoclonal antibody targeting human UAP56/DDX39B (Q13838), validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat. Recombinant immunogen (aa310-428). Ideal for mRNA export and spliceosome research.

[Learn More](#)

Parameter	Value
Product Name	UAP56/DDX39B Rabbit mAb
Remarks/Alias	BAT1; UAP56; D6S81E; UAP56/DDX39B
Species	Human
Gene ID (Human)	7919
Gene ID	7919
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-428 of human UAP56/DDX39B (Q13838).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13838
Protein Weight	49kDa
Shipping	Ice bag