


Anti-Uap56/Ddx39B Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13838

Item Number: CM0025422



Introduction

Rabbit polyclonal antibody against UAP56/DDX39B (BAT1), a 49 kDa spliceosome RNA helicase. Validated for WB, IHC-P, IF/ICC, and ELISA in human, mouse, and rat. Ideal for mRNA processing and export studies.

[Learn More](#)

Parameter	Value
Product Name	UAP56/DDX39B Rabbit pAb
Remarks/Alias	BAT1; UAP56; D6S81E; UAP56/DDX39B
Species	Human
GeneID (Human)	7919
GeneID	7919
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 319-428 of human UAP56/DDX39B (NP_004631.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13838
Protein Weight	49kDa
Shipping	Ice bag