


Anti-Uba1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P22314

Item Number: CM0017867



Introduction

High-quality anti-UBA1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat UBA1 (118 kDa). Ideal for ubiquitin-proteasome system studies.

[Learn More](#)

Parameter	Value
Product Name	UBA1 Rabbit pAb
Remarks/Alias	A1S9; A1ST; GXP1; UBE1; A1S9T; AMCX1; POC20; SMAX2; UBA1A; UBE1X; VEXAS; CFAP124; UBA1
Species	Human
Gene ID	7317
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-600 of human UBA1 (NP_695012.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22314
Protein Weight	118kDa
Shipping	Ice bag