


Anti-Uba1/7 Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P22314 / P41226

Item Number: CM0007828



Introduction

Recombinant rabbit monoclonal antibody targeting human UBA1 and UBA7 (UniProt P22314/P41226), validated for WB, IF-P, IHC-P, ELISA. Suitable for ubiquitin pathway research.

[Learn More](#)

Parameter	Value
Product Name	UBA1/7 Rabbit PolymAb®
Remarks/Alias	A1S9; A1ST; GXP1; UBE1; A1S9T; AMCX1; POC20; SMAX2; UBA1A; UBE1X; VEXAS; CFAP124
Species	Human
Gene ID (Human)	7317/7318
Gene ID	7317/7318
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22314/P41226
Protein Weight	118kDa/112kDa
Shipping	Ice bag