

Anti-Uba3 Rabbit Monoclonal Antibody For Wb, Elisa - Q8Tbc4

Item Number: CM0003660



Introduction

Rabbit monoclonal antibody against UBA3 (NEDD8-activating enzyme E1 catalytic subunit, 52 kDa). Validated for WB and ELISA; cross-reacts with human, mouse, and rat. Ideal for neddylation pathway and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	UBA3 Rabbit mAb
Remarks/Alias	NAE2, UBE1C, hUBA3, UBA3
Species	Human
GeneID (human)	9039
GeneID	9039
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-463 of human UBA3 (NP_003959.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TBC4
Protein Weight	52kDa
Shipping	Ice bag