

Anti-Uba5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Gzz9

Item Number: CM0014309



Introduction

Rabbit polyclonal antibody targeting UBA5 (Ubiquitin-like modifier-activating enzyme 5) for use in WB, IHC-P, and ELISA on human, mouse, and rat samples. Suitable for ufmylation studies.

[Learn More](#)

Parameter	Value
Product Name	UBA5 Rabbit pAb
Remarks/Alias	DEE44; EIEE44; SCAR24; THIFP1; UBE1DC1; UBA5
Species	Human
Gene ID (Human)	79876
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-404 of human UBA5 (NP_079094.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZZ9
Protein Weight	45kDa
Shipping	Ice bag