

Anti-Uba52 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62987

Item Number: CM0003467



Introduction

Rabbit monoclonal antibody against human UBA52 (P62987). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Detects ubiquitin-ribosomal fusion protein L40, suitable for ubiquitin pathway and translation research.

[Learn More](#)

Parameter	Value
Product Name	UBA52 Rabbit mAb
Remarks/Alias	L40; CEP52; RPL40; HUBCEP52; UBA52
Species	Human
GeneID (Human)	7311
GeneID	7311
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human UBA52 (P62987).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62987
Protein Weight	15kDa
Shipping	Ice bag