

Anti-Ubac2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nbm4

Item Number: CM0015365

Introduction

Rabbit polyclonal antibody targeting ubiquitin-associated domain-containing protein 2 (UBAC2). Validated for Western blot, immunofluorescence, and ELISA. Reacts with human, mouse, and rat. Involved in ER-phagy and Wnt signaling.

[Learn More](#)

Parameter	Value
Product Name	UBAC2 Rabbit pAb
Remarks/Alias	PHGDHL1; UBAC2
Species	Human
Gene ID (Human)	337867
Gene ID	337867
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-310 of human UBAC2 (NP_001137544.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NBM4
Protein Weight	39kDa
Shipping	Ice bag