

Anti-Ubac2 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Nbm4

Item Number: CM0023044

Introduction

High-quality Anti-UBAC2 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human and mouse UBAC2, an ER-phagy receptor and Wnt signaling regulator. Essential for autophagy and lymphocyte signaling studies.

[Learn More](#)

Parameter	Value
Product Name	UBAC2 Rabbit pAb
Remarks/Alias	PHGDHL1, UBAC2
Species	Human
GeneID (Human)	337867
GeneID	337867
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-310 of human UBAC2 (NP_001137544.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8NBM4
Protein Weight	39kDa
Shipping	Ice bag