

## Anti-Ubash3B Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Tf42

Item Number: CM0030893



### Introduction

Rabbit polyclonal antibody against UBASH3B (STS-1/TULA-2) for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying receptor tyrosine kinase signaling and ubiquitination pathways. 73kDa.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | UBASH3B Rabbit pAb                                                                                                 |
| Remarks/Alias    | p70, STS1, STS-1, TULA2, TULA-2, UBASH3B                                                                           |
| Species          | Human                                                                                                              |
| Gene ID          | 84959                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human UBASH3B (NP_116262.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q8TF42                                                                                                             |
| Protein Weight   | 73kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |