

Anti-Ube1L / Uba7 Monoclonal Antibody For Wb, Elisa - P41226

Item Number: CM0006532



Introduction

High-quality rabbit monoclonal antibody targeting human Ube1L/UBA7 for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying ISG15 conjugation and antiviral immunity. – 112kDa.

[Learn More](#)

Parameter	Value
Product Name	Ube1L / UBA7 Rabbit mAb
Remarks/Alias	D8, UBE2, UBE7, UBA1B, UBE1L, Ube1L / UBA7
Species	Human
Gene ID (Human)	7318
Gene ID	7318
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human Ube1L / UBA7 (NP_003326.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41226
Protein Weight	112kDa
Shipping	Ice bag