

Anti-Ube2A Polyclonal Antibody For Wb, Ihc-P, Elisa - P49459

Item Number: CM0028914



Introduction

Anti-UBE2A rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat UBE2A (RAD6A, 17kDa). Ideal for ubiquitin pathway, chromatin regulation, and mitophagy studies.

[Learn More](#)

Parameter	Value
Product Name	UBE2A Rabbit pAb
Remarks/Alias	UBC2; HHR6A; MRXSN; RAD6A; MRXS30; UBE2A
Species	Human
Gene ID (Human)	7319
Gene ID	7319
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-152 of human UBE2A (NP_003327.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49459
Protein Weight	17kDa
Shipping	Ice bag